

TRINITY FLUIDS, LLC

All Trinity Fluids, LLC brands are designed to promote positive results, if used properly. You can choose to not deviate from our instructions or you can experiment with it based on your level of experience and desired personal look and feel you desire.

Usage instructions for Cadaver Hydrating Fluid Eco Blend, CHF-EB

- Mix 1 **CHF-EB** into 28 H₂O (ratio 1:50). Blend continually for several minutes to ensure complete homogenous solution.
- Pour into spray or Yorker bottle.
- When ready to put cadaver away for extended periods apply topically **CHF-EB** all over areas in which you would like to maintain a soft and pliable condition.
- Pour a half gallon of mixed **CHF-EB** into a bucket. Use cheese cloth or non cotton cloth to soak in solution. Lay the soaked cloth over the cadaver extremities and torso to increase hydration.

Optional additives to solution.

- To KILL bacteria & mold growth **add 2 ounce of SANISOL 7** to every one gallon of mixed solution.