
TRINITY FLUIDS

Designed for Quality, Value and Professionals

All Trinity Fluids, LLC brands are designed to promote positive results, if used properly. You can choose to not deviate from our instructions or you can experiment with it based on your level of experience and desired personal look and feel you desire.

The following mixing instructions will give the embalmer a successful place to start.

PI 3N1 mixing instructions:

PI 3N1 Eliminates blood clots, increases lubricity in arteries, veins and capillaries for better circulation of arterial fluids. Also, removes heavy salts found in water softener systems.

PRE-INJECTION

- Bodies that have heavy clotting of blood, been in a cooler for more than 24 hours and or show signs of severe dehydration: 5 – 8 ounces per gallon.
- 1.0 to 1.25 gallons of total solution should be sufficient to remove heavy clotting.

CO-INJECTION

- Hard water issues use 5 - 9 ounces per 3 gallons of mixed Arterial Fluid.
- It is best to have your water tested to determine exact amount of PI 3N1 to be used. Ask your sales rep to help you with this task.
- Remove blockages, improve flow & lubricate the vascular system: 1 – 3 ounces per gallon
- Additional humectants required...use 2- 4 ounces per gallon of arterial solution.