
TRINITY FLUIDS

Designed for Quality, Value and Professionals

All Trinity Fluids, LLC brands are designed to promote positive results, if used properly. You can choose to not deviate from our instructions or you can experiment with it based on your level of experience and desired personal look and feel of the finished product.

The mixing instructions will give the embalmer a guarantee of a successful process.

SF Jaundice Fluid mixing instruction

- It is intended as a conditioner as to be able to use formaldehyde based fluids after the pre-injection is complete.
- 5 - 8 ounces per gallon of water. 3 gallons max. Do not use intermittent or closed drainage.
- Note- the tissue will become very firm.
- The fluid can be used as an arterial fluid. The fluid is rated as a 4 index. You may experience a softening of the tissues after 3 days. This does not mean the body is not preserved.
- Although in most cases additional embalming chemicals are not required, you can use additional AF series chemicals after completing the jaundice protocol. This should only be done if in your professional opinion desires a more aggressive treatment due to trauma or difficult case situations.
- **SPECIAL NOTE:** SF Jaundice can also be used in the final injection of any arterial fluid to give fast firmer action. Add no more than 8 ounces to the last gallon of arterial fluid. You may want to start out using less as to get the proper firming you desire.